

Medical_Image_Processing

This repository demonstrates the medical image processing according to the Task:1 'Image processing' for the selection process of Research assistant(RA) at Naamii.

Project Structure

```
.
├── data/           # Input CT scan data
├── results/
├── |-             # Output results and visualizations
├── src/           # Source code
│   ├── segment.py # Bone segmentation implementation
│   ├── contour_expansion.py # Contour expansion logic
│   ├── randomized_contour.py # Randomized contour adjustments
│   └── landmarks_detection_tibia.py # Landmark detection
├── utils/         # Utility functions
│   ├── nifti_file_handling.py # NIFTI file operations
│   ├── visualization.py      # Visualization
│   ├── animation.py          # Animation generation
│   └── voxels_operation.py    # Voxel manipulation
```

Usage

Run the main tasks using the following commands:

Task 1.1 - Bone Segmentation

```
python src/segment.py
```

Task 1.2 - Contour Expansion

```
python src/contour_expansion.py
```

Task 1.3 - Randomized Contour

```
python src/randomized_contour.py
```

Task 1.4 - Landmark Detection

```
python src/landmarks_detection_tibia.py
```

Results

The processing results are saved in the results/ directory:

- Task 1.1: Bone segmentation masks
- Task 1.2: Expanded contour masks
- Task 1.3: Randomized contour masks
- Task 1.4: Landmark coordinates and visualization
 - Submission: Final Submission of the task 1.4

Overview of task.

This project focuses on bone segmentation in 3D CT images volume. Task is to process a specific bone region as outlined in the task description and ultimately identify a key anatomical landmark.

The provided data is in NIfTI (Neuroimaging Informatics Technology Initiative) is a file format commonly used in medical imaging, particularly for MRI and CT scan data.

Data Analysis

- Volume Dimensions: Shape: (512, 512, 216) - This is a high-resolution volume with 512×512 pixels in each slice and 216 slices total
- Value Range: Range: [-3024.0, 1769.0] - This appears to be in Hounsfield Units (HU), which is standard for CT scans. Bone typically has values above 300-400 HU, while soft tissue is generally between -100 and 100 HU

From the MIP Images:

- The bones are clearly visible with good contrast -Both femur (thigh bone) and tibia (shin bone) are well-defined in the sagittal and coronal views
- The knee joint is clearly visible where they meet
- The circular view in the axial MIP indicates a good cross-sectional capture of the joint

From the Histogram

- There are distinct peaks in the histogram, indicating different tissue types
- The suggested threshold (15.00) seems to be between soft tissue and bone peaks
- There's a large peak at around -1000 HU (likely representing air)
- There appears to be another peak at around -3000 HU (could be outside the scan field)
- The bone intensities appear to be in the positive range, likely from about 200-1800 HU

Task 1.1 – Bone Segmentation

First task is to segment the femur and tibia regions from the provided CT image, as illustrated in Figure 1.B and save as nii.gz format. Apply only image processing techniques.

Based on Data analysis, here's how we should approach the bone segmentation:

Preprocessing

- Apply a Gaussian filter to reduce noise while preserving edges

Thresholding

- From the histogram, we can see that a threshold around 200-300 HU would likely separate bone from soft tissue
- After testing 200 HU is best **Morphological Operations**
- Fills small holes and smooths the mask using binary closing and custom operations (via `apply_morphological_operations` from `utils`).

Separation of Femur and Tibia

- After thresholding, we'll need to separate the two bones. This can be done using connected component analysis. Label connected components in the bone mask.
- Filter out small components based on a minimum size.
- There was only one component i.e. bone.
- Use 3D Bounding Box + Heuristic Z-Split: Extract bounding box of the largest component. Define a split plane at 40–50% of the bounding box height. Split the bone into upper (femur) and lower (tibia) based on this threshold. Assign the top half to femur and bottom to tibia.
- This approach is deterministic, simple, and more robust for knee CTs. Finally after experimentation the optimum is 50%.

Label connected components in the bone mask.

Post-processing

- Apply morphological operations (opening/closing) to remove small isolated regions and fill holes
- Use region properties to ensure we've correctly identified femur and tibia based on their shape and position

Output Generation

- Saves the segmentation results as NIfTI files (.nii.gz format) as required, creates three output files: Femur segmentation, Tibia segmentation and Combined segmentation with labels (1 for femur, 2 for tibia) **Visualization**
- Provides visualization functions to examine the segmentation results
- Shows slices from different anatomical planes
- Overlays segmentation on original data for easy assessment

Task 1.2 – Contour Expansion

- The second task requires to expand the segmented mask by 2 mm uniformly outward. 2mm should be parameter.

Key Features

Accurate Distance-Based Expansion

- Voxel Spacing Calculation: Extracts the real-world voxel dimensions from the NIfTI affine matrix
- Millimeter to Voxel Conversion: Converts the 2mm expansion distance to appropriate voxel units for each dimension
- Spherical Structuring Element: Creates a proper 3D spherical kernel for uniform expansion

Parameterized Expansion

- The expansion_mm parameter is configurable (default: 2.0mm) as required by the task
- Can easily be changed to any other value for different expansion distances

Robust Morphological Operations

- Uses binary_dilation with a custom spherical structuring element
- Handles anisotropic voxel spacing (different resolutions in x, y, z directions)
- Ensures truly uniform expansion in real-world coordinates

Overlap Handling

- Detects and handles overlapping regions between expanded femur and tibia masks
- Preserves original bone assignments in overlapping areas
- Provides detailed reporting of overlap statistics

Comprehensive Visualization

- Contour-based Visualization: Shows original contours in red and expanded contours in blue
- Multi-plane Views: Displays results in axial and coronal planes
- Before/After Comparison: Clear visualization of the expansion effect

Complete Pipeline Integration

- Saves results in NII.GZ format as required

Task 1.3 – Randomized Contour Adjustment

The third task involves randomizing the expanded segmented mask such that:

- The new mask lies between the original segmentation mask and the 2 mm expanded mask.
- The randomized mask must not exceed the 2 mm expansion limit.

- The contour must not shrink below the original segmentation contour.
- 2 mm and random value should be parameter.

Advanced Randomization Algorithm

- Distance-Based Probability: Uses distance transform to create realistic random variations
- Smooth Random Fields: Applies Gaussian smoothing to create coherent random patterns (not just noise)
- Probability Weighting: Combines distance-based and random factors for realistic results
- Multiple Random Seeds: Generates different randomized masks with different random seeds

Morphological Refinement

- Contour Smoothing: Applies morphological operations to smooth the randomized contours
- Overlap Resolution: Handles overlapping regions between femur and tibia intelligently
- Quality Control: Ensures the randomized masks maintain anatomical plausibility

Comprehensive Visualization

- Multi-Mask Comparison: Shows original (red), expanded (blue), and randomized (green/orange) contours
- Multiple Views: Displays results in both axial and coronal planes
- Clear Differentiation: Uses different colors for easy identification of different masks

Task 1.4 – Landmark Detection on Tibia

The final task focuses on the tibia. You are required to:

- Required to save the following segmentation masks in .nii.gz file:
 - Original Mask
 - 2 mm Expanded Mask
 - 4 mm Expanded Mask
 - Randomized Mask 1
 - Randomized Mask 2
- Identify the medial and lateral lowest points on the tibial surface, as demonstrated in Figure 1.D for all the above 5 masks. Please note that Figure 1.D is provided for illustration purposes only. You are not required to generate this image as part of your submission. Our team has tools for this visualization.
- Submit the coordinates of these points for evaluation we will compare your results with the ground truth.

Solution's Explanation

Coordinate Extraction

- Gets coordinates of all mask points using `np.where(mask > 0)`
- Converts to homogeneous coordinates for transformation

World Space Transformation

- Transforms voxel coordinates to real-world coordinates using the affine matrix
- This ensures measurements are in millimeters, not voxels

Lowest Points Detection

- Finds points within 5mm of the absolute lowest point
- Separates points into medial and lateral based on x-coordinate median
- Returns the lowest point from each region

Mask Generation

- Gets original masks through bone segmentation
- Creates 2mm and 4mm expanded masks
- Generates two randomized masks